

## Gibson Assembly

- Modified from Gibson *et al.* (2009)
- This assembly method is an isothermal, single-reaction method for assembling multiple overlapping DNA molecules. By coordinating the activity of a 5' exonuclease, a DNA polymerase and a DNA ligase two adjacent DNA fragments with complementary terminal sequence overlaps can be joined into a covalently sealed molecule, without the use of any restriction endonuclease.
- Preparation of DNA molecules for *in vitro* recombination.
- Generate the complementary sequence overlaps by PCR using the Phusion DNA polymerase. If necessary add 5 M Betaine to the reaction mix by reducing the amount of H<sub>2</sub>O to decrease the number of false PCR products.
- Identify the PCR products of interest by gel electrophoresis with known DNA standards.
- Extract the PCR products from the gel by cutting out the DNA fragments and clean them up by using a commercial gel clean-up kit.
- *in vitro* recombination:
- Assembly mixture:
  - 320  $\mu$ L 5x isothermal reaction buffer
  - 0.64  $\mu$ L of 10 U mL<sup>-1</sup> T5 exonuclease (for DNA molecules overlapping by greater than 150 bp add 3.2  $\mu$ L of 10 U mL<sup>-1</sup> T5 exonuclease)
  - 20  $\mu$ L of 2 U mL<sup>-1</sup> Phusion DNA polymerase
  - 160  $\mu$ L of 40 U mL<sup>-1</sup> taq DNA ligase
  - add ddH<sub>2</sub>O water up to a final volume of 1.2 mL
- Aliquot 15  $\mu$ L of the reagent-enzyme mix and store it at  $-20^{\circ}\text{C}$ .
- Thaw 15  $\mu$ L assembly mixture aliquot and keep it on ice until use.
- Add 5  $\mu$ L of the purified DNA molecules in equimolar amounts (between 10 and 100 ng of each DNA fragment).
- Incubate the resulting mixture at  $50^{\circ}\text{C}$  for 15 to 60 min, with 60 min being optimal.
- Transformation (via heat shock or via electroporation) without cleaning up the assembly product.

From: iGEM Bielefeld-CeBiTec