



Small-scale protein expression

- BL21 DE3 cells (New England BioLabs) were transformed with the expression plasmid. (use the transformation protocol).
- Overnight cultures of single colonies grown in LB media supplemented with maintenance antibiotics were diluted 1000-fold into fresh media (10 mL) with maintenance antibiotics and grown at 37 °C with shaking at 200 r.p.m. to OD600 ~0.4–0.6 before induction with 0.05-0.2% arabinose
- Keep one tube aside as uninduced control
- Cells were grown for a further 3 h at 37 °C with shaking at 230 r.p.m
- Cells were isolated by centrifugation at 5,000g for 10 min (Keep cell pellet)
- The resulting pellet was resuspended in 300 μ L 10mM TrisHCl pH 8+Lysozyme 0.1mg/ml PMSF 1mM mixed with 150mg of 100 μ m zirconia beads and lysed on a bead beater (3x 40s speed 6.0 with 1x 60s on ice in between incubated on ice) before centrifugation at 10,000g for 2 min (Keep supernatant)
- After separating the beads from the supernatant, another centrifugation follows 13,000g for 10 min (Keep both pellet AND supernatant)
- The supernatant was collected as the soluble fraction, and the resulting pellet was resuspended in an additional 300 μ L 10mM TrisHCl pH 8+Lysozyme 0.1mg/ml PMSF 1mM to obtain the insoluble fraction
- To 37.5 μ L of each fraction was added 12.5 μ L 4 \times Laemmli sample loading buffer (Bio-Rad) containing 2 mM dithiothreitol (DTT; Sigma-Aldrich)
- After vortexing, the fractions were incubated at 95 °C for 10 min
- 12 μ L of each soluble fraction and 6 μ L of each insoluble fraction was loaded into the wells of a cast gel (12%). 6 μ L of Protein Ladder (Thermo Scientific PageRuler Plus Prestained Protein Ladder) was used as a reference
- Samples were separated by electrophoresis at 180 V for 35 min in 1x SDS running buffer. Gels were stained with Coomassie Blue reagent for 30min (to overnight), then washed several times with H₂O (30min to overnight) before imaging with a gel dock